

Day1: NGS raw data → ASV

Marine eDNA Metabarcoding Analysis Using R

Jeju ARMS Case Study

Song, Chi-une

2026.08.31

Overview — workshop workflow

From a marine sample to ecological interpretation

01

Sampling & DNA Preparation

Collect and prepare the biological material




- 1 Marine sample collection
- 2 DNA extraction
- 3 COI barcode amplification

02

NGS Sequencing

Convert amplified DNA into sequencing reads



- 1 Library preparation
- 2 Illumina sequencing
- 3 Raw sequencing output

03

Bioinformatics Processing

Turn raw reads into reliable sequence variants



- 1 FASTQ quality control
- 2 Filtering & denoising
- 3 ASV inference with DADA2

04

Ecological Analysis

Connect sequence variants to community ecology



- 1 Taxonomic assignment
- 2 Alpha & beta diversity
- 3 Ecological interpretation

DAY 1 · Raw reads → ASVs

DAY 2 · ASVs → Ecology

WORKSHOP OVERVIEW

Jeju ARMS Case Study

Do different marine habitats support different ARMS-associated communities?

GJ

Gangjeong

Macroalgae-dominated habitat

Turf-forming algae

Depth: ~13 m


BM

Bomok

Coral-dominated habitat

Scleractinian corals

Depth: ~12 m



Study sites on the southern coast of Jeju Island

Lee et al. (2024)

CORE QUESTION

Different surrounding habitats → Different ARMS communities?

Visual survey vs. DNA metabarcoding

Two complementary approaches

Visual survey


- Morphological identification
- Observer-dependent
- Limited resolution for some taxa



Lee et al. (2024)

DNA metabarcoding

- Higher sensitivity: Detect cryptic organisms
- More reproducible
- Potentially higher taxonomic resolution



DNA metabarcoding can detect biodiversity overlooked by visual surveys.

01

DNA barcode & NGS sequencing

- What is a DNA barcode?
 - COI target region and primers
 - PCR product to sequencing library
 - Illumina paired-end sequencing
-

01

What is a DNA barcode?


Why do we target specific DNA regions?

Common barcode markers

- 16S rRNA — bacteria and archaea
- 18S rRNA — broad eukaryotic diversity
- COI — animals
- ITS — fungi

A useful DNA barcode should:

- be **amplifiable across many taxa**
- contain **enough sequence variation to distinguish taxa**
- have **well-represented reference database**



KEY POINT • A DNA barcode is a standardized genetic region used to distinguish and identify organisms.

COI target region and primers

Target marker: mitochondrial cytochrome c oxidase subunit I (COI)

Why COI?


- Widely used as the standard DNA barcode for animals
- Provides high taxonomic resolution for many animal groups
- Well represented in public reference databases

COI barcode region

- Mitochondrial protein-coding gene
- Full COI gene: approximately **1.5 kb**
- Standard animal barcode region: approximately **658 bp**

Primer pair used in this workshop


- **mlCOLintF (Forward)**
5' -GGWACWGGWTGAACWGTWTAYCCYCC-3'
- **jgHCO2198 (Reverse)**
5' -TAIACYTCIGGRTGICCRARAAYCA-3'



From PCR product to sequencing library

How do we prepare DNA fragments for Illumina sequencing?

- **PCR:** amplifies the target COI barcode region
- **Adapter:** allows DNA fragments to bind to the Illumina flow cell
- **Index:** identifies the sample of origin after sequencing




KEY POINT · PCR selects the target; adapters enable sequencing; indexes identify samples.


Illumina sequencing

How can millions of DNA fragments be sequenced at once?

- **Sanger sequencing:** sequences individual DNA templates separately
- **Illumina sequencing:** sequences millions of DNA fragments in parallel on a **flow cell**



Sanger sequencing




Illumina sequencing

KEY POINT • Illumina sequencing generates millions of reads in parallel.

Illumina sequencing

How does paired-end sequencing work?

- Both ends of the same DNA fragment are sequenced.
- **Read 1 (R1):** forward read
- **Read 2 (R2):** reverse read
- Overlapping regions can be merged into a single sequence.



Paired-end sequencing

KEY POINT • Paired-end sequencing reads the same DNA fragment from both ends, allowing overlapping reads to be merged.

02

BIOINFORMATICS PROCESSING

From raw sequencing reads to reliable ASVs

02

Raw sequencing output

What does the sequencing company actually give us?

Raw data statistics

Sample ID	Total read bases (bp)	Total reads	GC (%)	AT (%)	Q20 (%)	Q30 (%)
Bomok1	7,021,126	23,326	39.44	60.57	91.00	71.50
Bomok2	6,162,674	20,474	42.53	57.48	91.50	72.00
Bomok3	7,650,216	25,416	40.04	59.96	92.50	72.50
Bomok4	8,443,050	28,050	41.30	58.71	91.00	71.00
Gangjeong1	6,431,166	21,366	40.59	59.42	91.50	72.50
Gangjeong2	8,019,844	26,644	38.65	61.36	93.00	73.50
Gangjeong3	9,118,494	30,294	41.24	58.77	91.50	72.50
Gangjeong4	9,118,494	30,294	41.24	58.77	91.50	72.50

File contents

Download link	File size	md5sum
Bomok1_R1_001.fastq.gz	2.9 M	fb039b00b39bf46efe3b7db3bed9bf8b
Bomok1_R2_001.fastq.gz	3.2M	2cec4e1eede1a873a4fd24336a307033
Bomok2_R1_001.fastq.gz	2.6M	fad0e3e992afb749e4ffb3bd5dc11204
Bomok2_R2_001.fastq.gz	2.8M	62f0b05d654ce3907e1c5493aae669b6
Bomok3_R1_001.fastq.gz	3.2M	75615eb71d1309c74c8b2124b02a6c60
Bomok3_R2_001.fastq.gz	3.5M	5ffc12abf2dcb8bd7c4cf0a2e077a5ed
Bomok4_R1_001.fastq.gz	3.6M	9e4aeae7426ffee10d190da6c506266b
Bomok4_R2_001.fastq.gz	3.9M	80d7d19053393ea496abf3f0624e2f61
Gangjeong1_R1_001.fastq.gz	2.7M	f9f349e4f5a96202ecb49c88ad0cdf86
Gangjeong1_R2_001.fastq.gz	2.9M	73bcfe986778ec2a6043728ab7a624ed
...		

- Sequencing output

Sample A

→ **R1:** Sample_A_1.fastq.gz

→ **R2:** Sample_A_2.fastq.gz

- Sequencing statistics:

- read counts
- total bases
- GC/AT content

- Quality metrics: Q20 and Q30

KEY POINT · Raw sequencing output includes paired FASTQ files and basic sequencing-quality information.

What is inside a FASTQ file?

How are sequence and quality information stored?

- Each sequencing read is stored as a four-line FASTQ record.

FASTQ format

```
@SRR001666.1 071112_SLXA-EAS1_s_7:5:1:817:345 length=72      Field 1
GGGTGATGGCCGCTGCCGATGGCGTCAAATCCCACCAAGTTACCCTTAACAACCTTAAGGGTTTTCAAATAGA  Field 2
+SRR001666.1 071112_SLXA-EAS1_s_7:5:1:817:345 length=72      Field 3
|||||9IG9IC|||||ID|||||>|||||/                               Field 4
```

 Quality scores are encoded as ASCII characters.

Line	Contents
1	@ Read identifier
2	DNA sequence
3	+ Separator
4	Base-quality scores corresponding to each nucleotide in Line 2

KEY POINT • Every FASTQ read contains both sequence and quality information.


What is a Phred quality score?

How confident are we in each base call?

- Higher Q = higher confidence in the base call
 - Q20: ~99% base-call accuracy
 - Q30: ~99.9% base-call accuracy
- Modern Illumina FASTQ files commonly use Phred+33 encoding
- Read quality often declines toward the end of a read



Read1 (forward)



Read2 (reverse)



Why does this matter?

- Low-quality bases can create false sequence variants.

KEY POINT • Higher Phred scores mean more reliable nucleotide calls.

03

Practice

R/R studio setup
FASTQ quality control
ASV inference


03

R and RStudio

Why do we use them for metabarcoding analysis?




- **Free & open-source programming language and environment** for statistical computing
- Runs on **Windows, macOS, and Linux**
- Large ecosystem of packages for **bioinformatics, statistics, and ecology**
- Script-based analysis supports **reproducibility**



R Console



- An **integrated development environment (IDE)** for R
- Write and run R code directly from the source editor
- Organizes **scripts, files, plots, objects**, and **help** in one workspace
- **Projects** help organize analysis files and working directories



RStudio

03 PRACTICE

Installing R on Windows

What do we need before using RStudio?



1

DOWNLOAD

Download R for Windows → base
Click “[Download R-4.6.1 for Windows](#)”

* Note: macOS users should download the **macOS version**.

2

INSTALL

Run the Windows installer (.exe)
Keep the default installation options

3

CHECK

Open R and run: `R.version.string`
Confirm that R starts without errors

For this workshop: no PATH setup or Rtools installation is needed.

IMAGE PLACEHOLDER

Download and Install R

Precompiled binary distributions of the base system and contributed packages, **Windows and Mac** users most likely want one of these versions of R:

- [Download R for Linux](#) (Debian, Fedora/Redhat, Ubuntu)
- [Download R for macOS](#)
- [Download R for Windows](#)

R is part of many Linux distributions, you should check with your Linux package management system in addition to the link above.

<https://ftp.yz.yamagata-u.ac.jp/pub/cran/>

> `R.version.string` → Confirm the installed R version

```
> R.version.string  
[1] "R version 4.6.1 (2026-06-24 ucrt)"
```

KEY POINT • Install R first; RStudio will use the existing R installation.

03 PRACTICE

Installing RStudio Desktop

How do we get the interface for working with R?



1

DOWNLOAD

Go to the RStudio download page
Choose **“Download RStudio”** >> click the Download link

2

INSTALL

Run the Windows installer (.exe)
Keep the default installation options

3

LAUNCH

Open RStudio
Check that the Console starts with your installed R version

Important: R must be installed before RStudio. RStudio usually detects the installed R automatically.

IMAGE PLACEHOLDER

RStudio

The trusted IDE for R data scientists, now with AI that knows your data

The premier free, open-source development environment for R. Write code, explore data, and share results in one place, with built-in support for other languages R programmers use in their projects. Now with Posit AI for context-aware assistance.

[DOWNLOAD RSTUDIO](#) [TRY POSIT AI FREE](#)

[RStudio IDE | The Premier Free IDE for R](#)

RStudio IDE

Direct download links for RStudio Desktop 2026.08.1, the version this documentation describes:

Platform	Download
Windows 11	RStudio-2026.08.1-195.exe
Windows 11 (zip archive)	RStudio-2026.08.1-195.zip
macOS 14+	RStudio-2026.08.1-195.dmg
Ubuntu 22 / 24 / 26, Debian 13	rstudio-2026.08.1-195-amd64.deb
Ubuntu / Debian (tar.gz)	rstudio-2026.08.1-195-amd64-debian.tar.gz
RHEL 9 / 10, Fedora	rstudio-2026.08.1-195-x86_64.rpm

Troubleshooting

R should be installed before RStudio.

KEY POINT • RStudio is the workspace; it still requires R to perform the analysis.

03 PRACTICE

Meet RStudio

Where do we write code, see results, and manage files?

The screenshot displays the RStudio IDE with four numbered callouts indicating key components:

- 1 Source**: The Source pane on the left contains R code for locating paired FASTQ files, sorting them, and checking for mismatches. The code includes comments and function calls like `sort(list.files(...))` and `stop()`.
- 2 Console**: The Console pane at the bottom left shows the output of the executed code, including a list of file paths and a warning message about the 'package:stats' namespace.
- 3 Environment**: The Environment pane on the right lists the objects created during the analysis, such as `fnRs`, `FWD`, `max_ee`, `primer_dir`, `primer_mismatch`, `project_dir`, `raw_dir`, `result_dir`, `REV`, `sample_names`, `sample_names_F`, and `sample_names_R`.
- 4 Files / Plots**: The Plots pane on the right shows a quality score plot for the file `201908_ARMS_KJ_R2_001.fastq.gz`. The plot displays Quality Score (Y-axis, 0 to 40) versus Cycle (X-axis, 0 to 300). A red line indicates the quality score, and a red text label at the bottom left states "Reads: 10683".

Source

1 Write and save R scripts

Console

2 Run commands and view messages

Environment

3 Inspect objects created during analysis

Files / Plots

4 Manage files, plots, packages, and help

WORKSHOP TIP Run selected code with **Ctrl + Enter**

KEY POINT • Write code in the Source pane, run it in the Console, and inspect results in the other panes.

03 PRACTICE

RUN A QUICK CHECK

Load `dada2` and check that the `data/` folder is visible.

ARMS_Metabarcoding_Practice
유형: RStudio Project File

Download: ARMS_Metabarcoding_Practice

1 OPEN THE PROJECT

Open “ARMS_Metabarcoding_Practice.Rproj” before running the scripts

2 INSTALL THE PACKAGES

- scripts > open the “00_install_packages.R”
- Run all script code >> automatically installed

3 RUN A QUICK CHECK


Load `dada2` and confirm that the workshop data folders are visible

PACKAGE SETUP

```
install.packages(c("tidyverse", "vegan", "indicspecies"))
if (!requireNamespace("BiocManager", quietly = TRUE))
  install.packages("BiocManager")
BiocManager::install("dada2")
library(tidyverse)
library(vegan)
library(indicspecies)
library(dada2)

packageVersion("dada2")
```

IMAGE PLACEHOLDER



```
> packageVersion("dada2")
[1] '1.34.0'
```

KEY POINT • One RStudio Project + the same package set keeps the practical reproducible.

Why do raw reads need processing?

Why problems must we solve before inferring ASVs?

RAW SEQUENCING READS

Biological signal + technical noise

ATGCC G TAACTG	low Q
ATGCCGTAA C AG	error?
ATG T CGTAACTG	Chimeric?
C TATGCCGTAACTG	Adapters?
ATGCCGTAACTG	Identical?

Raw reads are observations — not yet biological sequence variants.

Repeated identical reads are useful abundance information; they can be collapsed computationally while preserving counts.

What must we handle?

- Low-quality bases**
 Trim or filter unreliable positions and reads.
- Sequencing / PCR errors**
 Model technical errors so they are not mistaken for real variants.
- Primer / adapter sequence**
 Remove non-biological sequence before downstream inference.
- Chimeric sequences**
 Identify PCR-derived combinations that do not represent true templates.
- Repeated identical reads**
 Dereplicate identical sequences for efficient computation while preserving read counts.

**QUALITY CONTROL
+ ERROR MODELING**

ASVs

KEY POINT

Raw reads contain both biological information and technical noise, so they must be processed before ASV inference.

How are ASVs different from OTUs?

Two different strategies for defining sequence units

OTU Cluster similar sequences

A predefined similarity threshold determines which sequences are grouped together.

Sequence A ATGCCGTA~~ACTG~~
Sequence B ATGCCGTAAC~~TA~~
Sequence C ATGCCGTAAC~~CG~~


} OTU 1

Example threshold: 97% similarity

The similarity threshold is chosen by the analyst.

ASV Infer sequence variants

DADA2 models sequencing errors and asks whether observed differences are consistent with error or with distinct variants.



ATGCCGTAAC × 124

ATGCCGTAAC ASV1

ATG~~T~~CGTAAC × 3


ATGTCGTAAC ASV2

No arbitrary similarity threshold is used to define ASVs.

03 PRACTICE

DADA2 workflow: from FASTQ to ASV table

How do raw sequencing reads become an analysis-ready sequence table?



TODAY'S FOCUS Understand why each step is performed — you do not need to memorize every function or parameter.

KEY POINT DADA2 progressively reduces technical noise while preserving biological sequence variation.

Practical 1 — Open the workshop project & settings

Let's make sure everyone starts from the same working environment.

1

OPEN THE PROJECT

Double-click ARMS_Metabarcoding_Practice.Rproj

2

CHECK THE FOLDERS

Confirm data/ scripts/ intermediate/ outputs/

3

OPEN THE SCRIPT

Open scripts/Day01_DADA2_ASV_inference.R

4

RUN THE FIRST CHECK


Run selected code with Ctrl + Enter

QUICK CHECK

Expected: R1/R2 FASTQ files

```
library(dada2)
list.files("data/fastq")
```

IMAGE PLACEHOLDER



When the setup is correct

```
> list.files("data/fastq")
[1] "Bomok1_R1_001.fastq.gz" "Bomok1_R2_001.fastq.gz"
[3] "Bomok2_R1_001.fastq.gz" "Bomok2_R2_001.fastq.gz"
[5] "Bomok3_R1_001.fastq.gz" "Bomok3_R2_001.fastq.gz"
[7] "Bomok4_R1_001.fastq.gz" "Bomok4_R2_001.fastq.gz"
[9] "Gangjeong1_R1_001.fastq.gz" "Gangjeong1_R2_001.fastq.gz"
[11] "Gangjeong2_R1_001.fastq.gz" "Gangjeong2_R2_001.fastq.gz"
[13] "Gangjeong3_R1_001.fastq.gz" "Gangjeong3_R2_001.fastq.gz"
[15] "Gangjeong4_R1_001.fastq.gz" "Gangjeong4_R2_001.fastq.gz"
```

KEY POINT

Open the RStudio Project first and confirm the input files before starting the analysis.

03 PRACTICE

Practical 2 — Inspect read quality

Where should we trim without losing paired-end overlap?

QUALITY

FILTER

ERROR

DENOISE

MERGE

CHIMERA

ASV TABLE

What are we looking for?

- Compare quality patterns between R1 and R2
- Find where quality declines near the read ends
- Choose trimming positions for the next step

R CODE



```
# 2. Inspect Read Quality (Raw) -----  
plotQualityProfile(fwd[1:2])  
plotQualityProfile(rev[1:2])
```

Do not trim by quality alone

Forward and reverse reads **must still overlap** after trimming.

Keep enough high-quality sequence for `mergePairs()`.

QUALITY PROFILE PLACEHOLDER



ASK THE CLASS

1. Which read loses quality earlier?
2. Where would you truncate R1 and R2?
3. Will enough overlap remain?

KEY POINT • Trimming must balance read quality and paired-end overlap.

Practical 3-5 — Filter and learn errors

How do we remove low-quality reads and model sequencing errors?

QUALITY

FILTER

ERROR

DENOISE

MERGE

CHIMERA

ASV TABLE

filterAndTrim()

Remove low-quality reads/bases using the trimming and filtering settings selected from the quality profiles.

learnErrors()

Learn the sequencing error pattern directly from the filtered data.

WORKSHOP CODE

```
# 3. Filter and Trim -----
primer_len <- c(26, 26) # COI FWD and REV primer lengths
trunc_len <- c(210, 180)

filt_f <- file.path(filt_path, paste0(sample_id, "_F.fastq.gz"))
filt_r <- file.path(filt_path, paste0(sample_id, "_R.fastq.gz"))

filtered <- filterAndTrim(
  fwd, filt_f, rev, filt_r, truncLen = trunc_len, trimLeft = primer_len,
  maxN = 0, maxEE = c(2, 4), truncQ = 2, rm.phix = TRUE, matchIDs = TRUE,
  compress = TRUE, multithread = FALSE)

# 4. Compare Before & After Trimming -----
plotQualityProfile(c(fwd[1], filt_f[1], rev[1], filt_r[1]))

# 5. Learn error rates -----
err_f <- learnErrors(filt_f, multithread = threads)
err_r <- learnErrors(filt_r, multithread = threads)

plotErrors(err_f, nominalQ = TRUE); plotErrors(err_r, nominalQ = TRUE)
```

Practical 3-5 — Filtering results and learned error model

What changed after filtering, and what error pattern did DADA2 learn?

QUALITY

FILTER

ERROR


DENOISE

MERGE

CHIMERA


ASV TABLE

filterAndTrim()



Low-quality tails are trimmed, and poor-quality reads are removed before denoising.

learnErrors()



DADA2 estimates how error rates change with quality score using the filtered data.

KEY POINT · Filtering removes unreliable reads, and error learning prepares the model used for denoising.

Practical 6 — Dereplication and denoising

How does DADA2 separate sequencing errors from real sequence variants?

QUALITY

FILTER

ERROR

DENOISE

MERGE

CHIMERA

ASV TABLE

derepFastq()

Collapse identical reads into unique sequences and record how many times each sequence occurs.

dada()

Use the learned error model to infer sequence variants separately from the forward and reverse reads.

WORKSHOP CODE

```
# 6. Dereplication and ASV inference (Denoise) -----  
# 1) Dereplication  
derepFs <- derepFastq(filt_f)  
derepRs <- derepFastq(filt_r)  
  
# 2) ASV inference  
asv_f <- dada(derepFs, err = err_f, multithread = FALSE)  
asv_r <- dada(derepRs, err = err_r, multithread = FALSE)
```

KEY POINT ·

Dereplication summarizes identical reads, and DADA2 uses the learned error model to infer real sequence variants.

Practical 7-9 — Merge reads, remove chimeras, and build the ASV table

How do denoised paired reads become the final ASV table?

QUALITY

FILTER

ERROR


DENOISE

MERGE**CHIMERA****ASV TABLE**

1 · Merge paired reads

Forward (R1)

Reverse (R2)



mergePairs()


Denoised R1 and R2 must overlap and agree in the shared region.

2 · Remove PCR chimeras

Parent A

Parent B

Chimera



removeBimeraDenovo()

Removes PCR artifacts formed by combining parts of different parent sequences.

FINAL PROCESSING

```
# 7. Merge paired reads -----
min_overlap <- 20

merged <- mergePairs(asv_f, derepFs, asv_r, derepRs,
                     minOverlap = min_overlap, maxMismatch = 0)

names(merged) <- sample_id
sequence_table <- makeSequenceTable(merged)

# 8. Remove chimeras -----
asv_matrix <- removeBimeraDenovo(
  sequence_table, method = "consensus", multithread = threads)

# 9. Export ASV table and representative sequences -----
abundance_order <- order(colSums(asv_matrix), decreasing = TRUE)
asv_matrix <- asv_matrix[, abundance_order, drop = FALSE]
...
```

KEY POINT · Denoised reads are merged, PCR chimeras are removed, and the remaining sequences form the final ASV table.

03 PRACTICE

What did DADA2 give us?

Output 1 — ASV abundance table

ASV × sample count matrix

Each value is the number of reads assigned to an ASV in a sample.

	Sample 1	Sample 2	Sample 3
ASV1	1250	21	0
ASV2	32	845	102
ASV3	0	15	520
...	...	...	...

ASV table

This table is the community matrix used for downstream ecological analysis.

How to read it

ROW

ASV

one inferred sequence variant

COLUMN

Sample

one biological sample

VALUE

Read count

abundance of that ASV
in the sample

Day 2 input: community abundance data

KEY POINT · The ASV table tells us how sequence variants are distributed across samples.

What did DADA2 give us?

Output 2 — Representative ASV sequences

Representative sequences (FASTA)

One nucleotide sequence is retained for each ASV.

```
>ASV1  
ATGGCACCTTTGACT...
```

```
>ASV2  
ATAGCTCGTACCAA...
```

```
>ASV3  
TTGCCGATGGTCAA...
```

Sequence identity only — no organism name yet

Do we know which organism ASV1 is?

Not yet.

DADA2 infers distinct sequence variants,
but it does not assign taxonomy.

**Sequence inference
≠ Taxonomic identification**

KEY POINT · Representative sequences are the starting point for BLAST and taxonomic assignment on Day 2.

Day 1 — What did we accomplish?

From ARMS samples to an ASV dataset

Today's workflow

Illumina

Paired-end sequencing

FASTQ

Raw reads + quality

QC + DADA2

Error-aware
sequence inference

ASV outputs

Table + sequences

Take-home messages

1

Raw reads contain technical errors

Quality control is required before treating sequence differences as biological variation.

2

DADA2 uses an error model

It learns sequencing-error patterns and infers ASVs rather than clustering at an arbitrary similarity threshold.

3

Day 1 ends with two key outputs

The ASV abundance table and representative ASV sequences become the starting point for taxonomy and ecology.

Overview — workshop workflow

From a marine sample to ecological interpretation

01

Sampling & DNA Preparation

Collect and prepare the biological material




- 1 Marine sample collection
- 2 DNA extraction
- 3 COI barcode amplification

02

NGS Sequencing

Convert amplified DNA into sequencing reads



- 1 Library preparation
- 2 Illumina sequencing
- 3 Raw sequencing output

03

Bioinformatics Processing

Turn raw reads into reliable sequence variants



- 1 FASTQ quality control
- 2 Filtering & denoising
- 3 ASV inference with DADA2

04

Ecological Analysis

Connect sequence variants to community ecology



- 1 Taxonomic assignment
- 2 Alpha & beta diversity
- 3 Ecological interpretation

DAY 1 · Raw reads → ASVs

DAY 2 · ASVs → Ecology

Day 1 answered “What sequence variants are present?” — Day 2 asks “What organisms are they, and how do communities differ?”